

# Storage under fume cupboards



### Application

The range of BECOME modules is designed, manufactured and certified in accordance with EN 14727, UNE-EN 16121 and UNE-EN 16122. It meets all their requirements, making it an ergonomic, safe product. Socket made of moisture-resistant material. It has a height-levelling system.

## Models



1. VG54/60 (PI, PD)



2. VG84 (P)

## Finishes


- Melamine.
- Fire resistant melamine.
- Water resistant melamine.
- Compact fronts.

## Colours


- White.
- Grey.

## Drawings

VG54/60



VG84



## Technical Characteristics

Reference	Model	Dimensions (mm)		
		Width	Depth	Height
VG54-PI	Left door	540	500	635
VG54-PD	Right door			
VG60-PI	Left door	600		
VG60-PD	Right door			
VG60-C3A	3 drawers			
VG84-P	Doors	840		