


Storage for acids under fume cupboards



Application

The range of BECOME modules is designed, manufactured and certified in accordance with EN-14727, UNE-EN 16121 and UNE-EN 16122. It meets all their requirements, making it an ergonomic, safe product. Cabinets for acid under fume cupboards have a removable shelf with polypropylene trays for a maximum load of 15 kg. The installation of a forced ventilation system by means of a polypropylene anti-corrosive extraction system is recommended.

Modelos



1. A27 VG54/60 (PI, PD)



2. A26 VG84 (P)

Finishes



- Melamine.
- Fire resistant melamine.
- Water resistant melamine.
- Compact fronts.

Colours



- White.
- Grey.

Drawings

A27 VG



A26 VG



Technical characteristics

Reference	Model	Dimensions (mm)		
		Width	Depth	Height
A27 VG54-PI	Left door	540	500	635
A27 VG54-PD	Right door			
A27 VG60-PI	Left door			
A27 VG60-PD	Right door			
A26 VG84-P	Doors	840		

Extraction diameter of 50 mm.