

Storage for acids in PP under fume cupboards



The range of BECOME modules is designed, manufactured and certified in accordance with EN 14727 standard. It meets all their requirements, making it an ergonomic, safe product. Made of solid panels and polypropylene components.. Removable storage shelf with polypropylene trays with a maximum load of 30 kg. Capacity to retain fluids in the event of spillages. 5 litres. The installation of a forced ventilation system by means of a polypropylene anti-corrosive extraction system is recommended.

Models



1. A27 PP VG54/60 (PI, PD)



2. A27 PP VG84 (P)

Finishes

- Polypropylene.

Colours



- Grey.

Accessories



- Extraction equipment.
- Filtration - ventilation box.

Drawings

A27 VG PP



A26 VG PP



Technical data

Reference	Model	Dimensions (mm)		
		Width	Depth	Height
PP A27 VG54-PI	Left door	540	500	635
PP A27 VG54-PD	Right door			
PP A27 VG60-PI	Left door	600		
PP A27 VG60-PD	Right door			
PP A26 VG84-P	Doors	840		

Extraction diameter of 75 mm.