

Storage for solvents under fume cupboards



Application

The range of safety cabinets is designed, manufactured and certified in accordance with EN-14727. It meets all their requirements, making it an ergonomic, safe product. Type 90 classification in accordance with EN 14470- 1. Metal body made of steel plate with a plastic paint powder coating. Insulating filler composed of several layers of mineral materials Intumescent gaskets for sealing all of the gaps and spaces between the door and the body, which expand in the event of fire and prevent the entry of heat into the cabinet. Series earth connection on the rear wall of the cabinet.

Models



1. S 30A pull-out drawer

2. S 31/33A with two

3. S 32A with three

Finishes


- Metal.

Colours


- Grey.

Drawings


S 30




S 32



S 33



S-31



Technical data

Reference	Model	Dimensions (mm)		
		Width	Depth	Height
S 30	Pull-out drawer	600	595	635
S 31	2 doors	1.100		
S 32	3 doors	1.400		
S 33	2 doors	888		