

RANGE 3030.T**EN 14470-1 FM 6050**

EN 14470-1 & FM⁽¹⁾ SAFETY CABINETS

TYPE 30 MINUTES - For flammable products with ergonomic retention drawers



*3034T5



*3035T5

CONFORMITY

- Tested and certified type 30 minutes according to EN 14470-1 (04/2004) and EN 1363-1 (06/2000).
- Tested and certified 10 minutes according to FM 6050 (Factory Mutual international).
- Normalized labelling according to ISO 3864 and ISO 7010, NF X08.003, Directive 92/58 / EEC, EN 14470-1 (04-2004) and FM specifications.

ADVANTAGES

- 30-minutes fire resistance according to ISO 834 fire heating curve.
- Pedestal area under driveable by a lorklift by removing panel
- Self-closing cabinets, with predominant wings and doors in open position thanks to hot melt fuses starting at 47 °C.
- Many models available.
- Suitable for all types of chemicals commonly used in laboratory and industry.
- Separation by families of hazardous products in accordance with the Public Health Code and the recommendations of the I.N.R.S ..
- Self-closing retention drawers (1) in case of fire

PASSIVE SECURITY

- Construction in insulating panels offering a high chemical resistance to aggressive products.
- Light gray frame RAL 7035 with yellow plated fields RAL 1023.
- Signaling with standardized symbols.

ACTIVE SECURITY

- Ventilation air outlet and inlet Ø 100 mm for possible connection.
- Self-closing doors with predominant wings.
- Self-closing device of the doors in case of fire.
- 3-points lock.
- Anti spark lock closing.
- Adjustable leveling feet.
- In case of fire, a thermo-fuse automatically seals off the ventilation system.
- Thermo-expanding seals of the door(s) allow keeping the internal body fireproof.

OPTIONS

- Additional drawers



Reference	Description	Type	Number of doors	External dimensions H x W x D (mm)	internal dimensions H x W x D (mm)	Weight (kg)
3034T4	Cabinet 1 door 4 drawers	High	1	1980 x 605 x 550	1580 x 495 x 420	141
3034T5	Cabinet 1 door 5 drawers	High	1	1980 x 605 x 550	1580 x 495 x 420	147
3034T6	Cabinet 1 door 6 drawers	High	1	1980 x 605 x 550	1580 x 495 x 420	153
3035T4	Cabinet 2 doors 4 drawers	High	2	1980 x 1115 x 550	1580 x 1050 x 420	258
3035T5	Cabinet 2 doors 5 drawers	High	2	1980 x 1115 x 550	1580 x 1050 x 420	270
3035T6	Cabinet 2 doors 6 drawers	High	2	1980 x 1115 x 550	1580 x 1050 x 420	287
T148	Additional drawer for 1 door cabinet ⁽¹⁾	High	-	80 x 420 x 380	-	6
T235	Additional drawer for 2 doors cabinet ⁽¹⁾	High	-	80 x 930 x 380	-	12

⁽¹⁾ To be ordered during manufacture. Maximum quantity per cabinet : 6

Flammable cabinets

Corrosive cabinets

Toxic cabinets and multirisk

Pesticide range

Fume hoods, filtering cabinets - Ventilation

Containment and Cans

File cabinets Anti-fire equipments

Showers and first-aid equipments