

RANGE 790+T

EN 14470-1 EN 1363-1 FM 6050



EN 14470-1 SAFETY CABINETS TYPE 105 MIN FOR FLAMMABLE PRODUCTS WITH SELF ERGONOMIC RETENTION DRAWERS



794+T5 ▲
Self retracting drawer with automatic closing in case of fire⁽²⁾



▲ 795+T6



792T ▲

P/N	Description	Type	Number of doors	Outside Dimensions H x W x D (mm)	Inside Dimensions H x W x D (mm)	Weight (kg)
791T	Cabinet 1 drawer	under-bench	1	670 x 635 x 615 ⁽¹⁾	450 x 480 x 370	155
792T	Cabinet 1 drawer	under-bench	1	670 x 1137 x 615 ⁽¹⁾	450 x 950 x 370	240
794+T4	Cabinet 1 door, 4 drawers	tall	1	1950 x 635 x 620	1620 x 490 x 410	301
794+T5	Cabinet 1 door, 5 drawers	tall	1	1950 x 635 x 620	1620 x 490 x 410	307
794+T6	Cabinet 1 door, 6 drawers	tall	1	1950 x 635 x 620	1620 x 490 x 410	313
795+T4	Cabinet 2 doors, 4 drawers	tall	2	1950 x 1137 x 620	1620 x 1000 x 410	483
795+T5	Cabinet 2 doors, 5 drawers	tall	2	1950 x 1137 x 620	1620 x 1000 x 410	495
795+T6	Cabinet 2 doors, 6 drawers	tall	2	1950 x 1137 x 620	1620 x 1000 x 410	507
T48	Additional drawer 1 door ⁽²⁾	tall	2	80 x 420 x 380	-	6
T35	Additional drawer 2 doors ⁽²⁾	tall	2	80 x 930 x 380	-	12

(1) Depth with the connecting sleeve at the back. Depth without connecting sleeve : 565 mm
(2) To be ordered during manufacture. Maximum quantity per cabinet: 6